

SALTER LABS®



Salter-Style®

Original Cannulas

The Original Salter-Style® Cannula is for oxygen flows up to 6 LPM. The clear, lightweight adult cannula is made with the Salter-Style® manufacturing process, producing a soft, anatomically designed, one piece, seamless cannula that is comfortable to wear.



Features

- ☐ Anatomically curved to fit upper lip
- ☐ Curved, tapered nasal prongs
- ☐ Over-the-ear style headset
- ☐ Three-channel safety tubing
- ☐ Soft, lightweight, clear material

Benefits

- ☐ Eliminates flap or ridge which can cause irritation
- ☐ Improved anatomical fit for increased comfort
- ☐ Secure positioning
- ☐ Uninterrupted flow of oxygen
- ☐ Better comfort for long-term use

Wrap around, "Anti-Roll" Feature

Works with the anatomy rather than against it...A Salter Innovation.




slc-03

Salter-Style® Original Cannula P/N 1600

Clear, lightweight cannula with unique, one piece anatomical design eliminates the need for flap, which can cause irritation. The result is maximum patient comfort. Over-the-ear style headset ensures secure positioning, allowing freedom of movement. Provided with three-channel safety tubing or smooth bore tubing in lengths from 4 to 100 feet, other special lengths available.

One piece,
anatomical
design




Three-channel
safety tubing

Conventional Cannula P/N 1056

Conventional cannula with non-flared tips. Made of soft material. Adjustable tubing fits comfortably over the ears to hold the cannula firmly in place. Provided with standard smooth bore tubing in 7' or 16' lengths. Special lengths available. See price list for details.

Non-flared
tips

Smooth
bore or
three-
channel
safety
tubing



Ordering Information

Salter-Style® Original Cannulas*

	Units (per case)	Safety Channel	Smooth Bore
Clear cannula with 7' supply tube	50	1600 ^x	1650
Clear cannula without supply tube	50	1606 ^x	
Clear cannula with 7' supply tube and female thread grip connector	50	1699	

Special Length Salter Cannulas*

Adult cannula with 1' supply tube	50	1600-1	
Adult cannula with 4' supply tube (ideal for strollers and portable units)	50	1600-4 ^x	1650-4
Adult cannula with 9' supply tube	50	1600-9	1650-9
Adult cannula with 10' supply tube	50	1600-10	1650-10
Adult cannula with 12' supply tube	50	1600-12	1650-12
Adult cannula with 13' supply tube	50	1600-13	1650-13
Adult cannula with 14' supply tube	50	1600-14 ^x	1650-14
Adult cannula with 15' supply tube	50	1600-15 ^x	1650-15
Adult cannula with 16' supply tube	50	1600-16	1650-16
Adult cannula with 20' supply tube	25	1600-20	1650-20
Adult cannula with 21' supply tube	25	1600-21 ^x	1650-21
Adult cannula with 24' supply tube	25	1600-24	1650-24
Adult cannula with 25' supply tube	25	1600-25 ^x	1650-25
Adult cannula with 30' supply tube	20	1600-30	1650-30
Adult cannula with 35' supply tube	20	1600-35	1650-35
Adult cannula with 40' supply tube	20	1600-40	1650-40
Adult cannula with 50' supply tube	20	1600-50	1650-50
Adult cannula with 75' supply tube	15	1600-75	1650-75
Adult cannula with 100' supply tube	10	1600-100	1650-100

Special lengths available on request. Minimum order is one (1) case.

Conventional Style Cannulas

Conventional cannula with 7' supply tube	50	1053 ^{x*}	1056 [*]
Conventional cannula without tube	50	1057 ^x	
Conventional cannula with 16' supply tube	50	1063 [*]	1068 [*]

^x This cannula is available with closed cell, over-the-ear, foam cushions, by adding TLC to the part number. (see catalog page SLC-00)

SALTER LABS®

100 W. Sycamore Road - Arvin, CA 93203
Phone: 661-854-3166
1-800-421-0024
1-800-235-4203
Fax: 661-854-3850
Toll Free Fax: 1-800-628-4690
Fax: 661-854-3850
E-mail: salterlabs@us.salterlabs.com
www.salterlabs.com

CE
0482