


SALTER LABS®



Salter Labs®

Elongated Medium Concentration Oxygen Masks

Two Salter advantages... crush-resistant tubing plus the ease and security of an over-the-ear style.



Features

- ☐ Transparent mask with soft anatomical form
- ☐ Adult mask in over-the-ear style . A Salter exclusive innovation!
- ☐ All masks come standard with three-channel safety tubing
- ☐ Wide, adjustable elastic strap also available
- ☐ Elastic strap ends sealed

Benefits


- ☐ Better fit for long-term users and edentulous patients
- ☐ Fits comfortably over the ears. Eliminates need to move patient's head to apply. Useful for head, neck and spinal case patients
- ☐ Crush-resistant safety tubing reduces potential for occlusion of oxygen supply
- ☐ Fits comfortably below the ears to hold mask snugly in place
- ☐ Precludes band being pulled through mask



slc-51

Elongated Medium Concentration Oxygen Mask

Transparent mask with soft anatomical form ideal for long-term users and edentulous patients. Over the ear style fits comfortably to hold the mask snugly in place. Exceptionally comfortable for even hard-to-fit faces. Three-channel crush resistant tubing standard.



P/N 8005



Ordering Information

Medium Concentration Oxygen Masks

	Units (per case)	Elastic-Strap Style	Over-The-Ear Style
Adult elongated medium concentration oxygen mask with 7' safety tube	50	8110	8005
Adult Elongated medium concentration oxygen mask	50	8115	8010
Adult extra-large medium concentration oxygen mask with 7' safety tube	50	1102	

Pediatric medium concentration oxygen mask with 7'safety tube	50	1122
Pediatric medium concentration oxygen mask without tube	50	1123
Infant medium concentration mask with 7' safety tube	10 or 50	1114
Infant medium concentration mask without safety tube	10 or 50	1115



100 W. Sycamore Road - Arvin, CA 93203
 Phone: 661-854-3166
 1-800-421-0024
 1-800-235-4203
 Fax: 661-854-3850
 Toll Free Fax: 1-800-628-4690
 Fax: 661-854-3850
 E-mail: salterlabs@us.salterlabs.com
www.salterlabs.com



0482